

# Immediate Biomarker Profiling of Actionable NSCLC Panel from a Digitized H&E Biopsy Image

Nurit Paz-Yaacov<sup>1</sup>, Efrat Ofek<sup>2</sup>, Razan Haj<sup>2</sup>, Yossef Molchanov<sup>2</sup>, Rinat Yacobi<sup>2</sup>, Chen Mayer<sup>2</sup>, Camila Avivi<sup>2</sup>, Jonathan Zalach<sup>1</sup>, Inbal Gazy<sup>1</sup>, Nir Peled<sup>1,3</sup>, Jair Bar<sup>4</sup>, Dov Hershkoviz<sup>5,6</sup> and Iris Barshack<sup>2,6</sup>

<sup>1</sup> Imagene AI LTD., <sup>2</sup> Institute of Pathology, Sheba Medical Center, <sup>3</sup> Department of Oncology, Shaare Zedek Medical Center, <sup>4</sup> Institute of Oncology, Sheba Medical Center, <sup>5</sup> Pathology Department, Tel Aviv Sourasky Medical Center, <sup>6</sup> Sackler Faculty of Medicine, Tel Aviv University

## Background

- Optimal non-small cell lung cancer (NSCLC) patient management rely on comprehensive molecular testing to identify driver alterations that guide therapeutic decisions.
- Limited accessibility, long turnaround times and inadequate tissue for molecular testing all result in undertesting and suboptimal clinical management of lung cancer patients.
- An AI solution that uses the diagnostic H&E-stained slide to detect genomic alterations offers an accessible solution that can overcome many of the aforementioned barriers to molecular testing.

## Aim

To evaluate the accuracy of an AI-based solution that uses the diagnostic H&E-stained slide to detect the NSCLC actionable genomic alterations compared to clinically diagnostic tests.

## Method

- Training and testing was performed on pan-cancer FFPE H&E digital whole slide images (WSIs) from both public and commercial databases as well as medical centers using self-supervised learning (SSL) with dynamic data augmentation and multiple instance learning (MIL) algorithms.
- The model performance was evaluated on an independent lung cancer WSI dataset (including samples from multiple sites, led by Sheba Medical Center) using a 5-fold cross-validation scheme, and comparing the inference results to the official reported results of the reference methods (NGS, FISH, etc.) for each biomarker.
- For every sample, each gene was assigned either a "positive" or a "negative" status for high confidence biomarker scores, or "inconclusive" status for low confidence biomarker scores.

## Data

- 1,915 retrospective samples** were used for the evaluation of the model:
  - Multiple medical centers and reference labs, from the USA and Israel, led by the Pathology Institute at Sheba Medical Center.
  - Various digital slide scanners, including Philips and Leica devices.
  - Needle biopsies as well as surgical resection specimens.
  - Specimens from both primary and metastatic sites.
- All collected samples from NSCLC patients with an available WSI and an official result were included in the study. Each case was represented by a single slide.

## Results

AI-based NSCLC actionable biomarker panel performance:

Biomarker	N cohort*	N pos	N neg	AUC	Sensitivity	Specificity	Accuracy	% of cohort	Prevalence
<i>ALK fusion</i>	1513	88	1425	0.979	0.974	0.962	0.962	53	0.05
<i>BRAF V600E</i>	1748	43	1705	0.983	0.941	0.973	0.973	65	0.02
<i>EGFR mut</i>	1870	318	1552	0.973	0.901	0.917	0.914	44	0.15
<i>ERBB2 mut</i>	1642	46	1596	0.990	1.000	0.973	0.973	43	0.03
<i>MET ex14 skipping</i>	1330	67	1263	0.946	0.952	0.921	0.923	25	0.04
<i>NTRK1-3 fusion</i>	1082	25	1057	0.852	0.824	0.944	0.942	94	0.01
<i>RET fusion</i>	1308	43	1265	0.941	0.923	0.920	0.920	55	0.02
<i>ROS1 fusion</i>	1304	46	1258	0.984	1.000	0.970	0.970	51	0.02
<b>Average</b>				<b>0.956</b>	<b>0.939</b>	<b>0.947</b>	<b>0.947</b>	<b>54</b>	

\* Number of samples with official results for each biomarker.


*KRAS* G12C biomarker was not included in the panel yet as it's performance values were: 91.7% specificity, 82.4% sensitivity, 91.2% accuracy and an AUC of 0.842 for ~0.25 of the cohort's patients.

NSCLC inference performance, in the high-confidence cases, demonstrated an **average of:**

AUC	Sensitivity	Specificity	Accuracy
<b>0.956</b>	<b>93.9%</b>	<b>94.7%</b>	<b>94.7%</b>

## Discussion

- The NSCLC panel performances demonstrated **generalization** in detecting genomic alterations directly from diagnostic H&E-stained WSIs with **high accuracies**, comparable to those of diagnostic tests.
- We demonstrate that the NSCLC actionable alterations are **identifiable using image-based AI algorithms**. Imagene is currently researching the boundaries of AI solutions capabilities such as the ability to distinguish between different variants within a specific gene.
- Further studies evaluating the practical benefits of **incorporating NSCLC panel into the clinical workflow** are currently underway. Such studies will assess how a rapid, robust, scalable AI biomarker profiling panel can optimize and improve the diagnostic workflow and enhance accessibility for molecular testing.



Schematic representation in the diagnostic workflow